
TECHNICAL INFORMATION

Indented Writing Restoration Kit Catalog No. IRS1001

INTRODUCTION

The IRS1001 is primarily recommended for developing invisible writing, and also for rendering paper indentations and abrasions visible. Whenever writing is made on a leaf or sheet of a pad or book, not only is the top sheet affected, but the under sheets also receive impressions from which it is possible to determine the text of the original writing. These impressions are called “indented writing” and it is possible to develop these latent impressions with the reagents furnished with this kit. Individual application vials for the reagents provide unlimited storage life.



100 HUNTER PLACE, YOUNGSVILLE, NC 27596 USA

Ph: (919) 554-2244, (800) 356-7311 • Fax: (919) 554-2266, (800) 899-8181 • Web: www.sirchie.com • Email: sirchieinfo@sirchie.com