



# Narcotics Investigation

As the opioid epidemic grows and becomes more dangerous to law enforcement and first responders, our goal is to make you as safe as possible. We've reviewed the Drug Enforcement Administration (DEA) and the Centers for Disease Control (CDC) guidelines for handling fentanyl, carfentanil and acetyl-fentanyl.

Building from guidelines offered by the Drug Enforcement Administration (DEA), the Centers for Disease Control (CDC), and the manufacturers of pharmaceutical fentanyl, we offer solutions under our NARKsafe™ product line.

NARKsafe™ solutions are designed to help keep you and your team safe while handling, packaging, transporting, testing, and storing potential fentanyl, carfentanil, acetyl-fentanyl, and other dangerous substances.

For more information about these guidelines, please see:

- Fentanyl: A Briefing Guide for First Responders, DEA [https://www.dea.gov/druginfo/Fentanyl\\_BriefingGuideforFirstResponders\\_June2017.pdf](https://www.dea.gov/druginfo/Fentanyl_BriefingGuideforFirstResponders_June2017.pdf)
- Fentanyl: Preventing Occupational Exposure to Emergency Responders, CDC <https://www.cdc.gov/niosh/topics/fentanyl/default.html>

NARKsafe™



Fentanyl Facts:

- Fentanyl is a synthetic opioid developed as pain management treatment
- 50 to 100 times more potent than morphine and 25 to 40 times more potent than heroin
- Appearance: usually white powder similar to methamphetamine and cocaine
- Fentanyl is a Schedule II substance (Heroin is Schedule I)
- Attractive to drug abusers due to its intense high
- Street names: China Girl, China White, Goodfella, Apache, TNT, Tango & Cash, Dance Fever, Jackpot
- Trade Names: Actiq®, Fentora™, Duragesic®, Sublimaze®
- Clandestinely produced and either cut into heroin or mixed with other adulterants and sold as heroin
- Signs of opioid overdose: pinpoint pupils, respiratory depression and unconsciousness
- Antidote: Narcan (Naloxone) availability: <https://www.bjatraining.org/tools/naloxone/Naloxone-Background>

Carfentanil Facts:

- 100 times more powerful than fentanyl
- Medically used to sedate large animals such as buffalo and elephants
- Two milligrams of carfentanil can sedate a 2,000 pound African elephant



Often Bought Together:  
Indicates items often purchased in conjunction with other items.

NARK II Presumptive Drug Testing Reagent Pouches

Our NARK drug identification pouches can presumptively identify several families of substances that could be abused drugs. Each of the tests is comprised of one or more chemical reagents. When a color or series of colors occurs within a specific testing sequence, a positive identification may be presumed.



Hemp/CBD Detection Test (NARK20034)

This latest addition to our NARK II Product Line (NARK20034) is a new screening test that enables you to presumptively determine if an unknown plant or oil substance has suspect low levels of THC and is a legal hemp/CBD product, or if it has suspect high levels of THC, and is marijuana. If the substance tests positive for marijuana, use our Duquenois-Levine Reagent (NARK2005) or our KN Reagent (NARK2020) to confirm.



This new screening test saves you the hassle, the embarrassment, and the time associated with arresting/booking someone for possession of legal Hemp or CBD products.

No.	Description	Presumptively tests for
NARK2001	Marquis Reagent, 10/box	Opium Alkaloids (Heroin/Morphine), MDMA (Ecstasy), Amphetamines/Meth
NARK2002	Nitric Acid Reagent, 10/box	Differentiate Heroin/Morphine
NARK2003	Dille-Koppanyi Reagent, 10box	Barbiturates
NARK2004	Ehrlich's Reagent, 10/box	LSD
NARK2005	Duquenois-Levine Reagent, 10/box	Marijuana, Hash, Hash Oil, THC Concentrates
NARK2007	Scott Reagent (Modified), 10/box	Cocaine Salts (HCl), Cocaine Base (Crack)
NARK2008	Methadone Reagent, 10/box	Methadone
NARK2009	PCP/Methaqualone Reagent, 10/box	PCP, Methaqualone
NARK20010	Special Opiates Reagent, 10/box	Codeine, Heroin, Morphine, Oxycodone, Buprenorphine
NARK20011	Mecke's Modified Reagent, 10/box	Heroin (White, Brown, Black Tar)
NARK20012	Frohdes Reagent, 10/box	Talwin® (Pentazocine), DMT (Dimethyltryptamine)
NARK20013	Ephedrine Reagent, 10/box	Ephedrine, Psuedoephedrine
NARK20014	Valium Reagent, 10/box	Valium, Rohypnol, Ketamine
NARK20015	Methamphetamine/MDMA Reagent, 10/box	Methamphetamine, MDMA (Ecstasy)
NARK20019	Mayer's Reagent, 10/box	General Screening, Narcotic Alkaloids
NARK20020	KN Reagent (Fast Blue B Salts), 10/box	Marijuana, Hashish, Hash Oil, THC Concentrates (DAB, WAX, BHO)
NARK20021	GHB Reagent, 5/box	GHB (Gamma Hydroxybutyrate)
NARK20022	Mandelin Reagent, 10/box	General screening for Methadone, Amphetamine-Type Compounds
NARK20023	Synthetic Cannabinoids Reagent, 5/box	Synthetic Cannabinoids (indole formulations)
NARK20024	MDPV Reagent, 10/box	MDPV, Butylone, Ethylone (bk-MDEA), Methylone, Naphyrone, Pentylone
NARK20025	Mephedrone Reagent, 10/box	Mephedrone (Synthetic Cathinones)
NARK20026	a-PVP Reagent, 10/box	a-PVP (Flakka)
NARK20029	2C Reagent, 10/box	2C-B, 2C-B-NBOMe, 2C-C, 2C-C-NBOMe, 2C-E, 2C-E-NBOMe, 2C-I, 2C-I-NBOMe, 2C-N, 2C-T-7
NARK20030	Psilocybin/Psilocin Reagent, 10/box	Psilocybin, Psilocin (Mushrooms)
NARK20031	Liebermann Reagent, 10/box	Methamphetamine, Morphine
NARK20032	Mollies Reagent, 10/box	Mollies
NARK20033	Fentanyl Reagent, 10/box	Fentanyl, Acetyl-fentanyl, Heroin
NARK20034	Hemp/CBD Detection Test, 10/box	Hemp, Cannabinoids/CBD
No.	Description	
NARK2006	Acid Neutralizer Liquid, 1 oz bottle	
NARK200N	Acid Neutralizer Powder, 0.5 oz. pouch, 10 ea.	



**Pricing Information:**  
For current pricing visit  
[www.sirchie.com](http://www.sirchie.com), or  
call 800.356.7311

**Often Bought Together:**  
Indicates items often  
purchased in conjunction  
with other items.

NARK II Notebook and Patrol Kits

Carry the most commonly used tests with you in a protected case, keeping the glass ampoules in the pouches safe. We designed our Patrol Kit with the reagents most often encountered by patrol officers.



NARK12PAK Contents:

- 2- NARK2001-1 Marquis Reagent (opium alkaloids)
- 2- NARK2005-1 Duquenois-Levine Reagent (marijuana)
- 4- NARK2007-1 Scott Reagent Modified (cocaine HCl/crack)
- 2- NARK20011-1 Mecke's Reagent Modified (all heroin)
- 2- NARK20015-1 Sodium Nitroprusside Reagent (methamphetamine/MDMA)
- 1- NARK12PAKC Copolymer Case, Black, 12.5" x 10.25" x 1.25"

NARK24PAK Contents:

- 4- NARK2001-1 Marquis Reagent (opium alkaloids)
- 6- NARK2005-1 Duquenois-Levine Reagent (marijuana)
- 6- NARK2007-1 Scott Reagent Modified (cocaine HCl/crack)
- 4- NARK20011-1 Mecke's Reagent Modified (all heroin)
- 4- NARK20015-1 Sodium Nitroprusside Reagent (methamphetamine/MDMA)
- 1- NARK24PAKC Copolymer Case, Black, 12.5" x 10.25" x 1.25"

NARK200VC Contents:

- 2- NARK2005-1 Duquenois-Levine Test for Marijuana, Hashish, Hash-Oil, THC
- 4- NARK2007-1 Scott Reagent (modified) Test for Cocaine Salts, Cocaine Base
- 2- NARK20015-1 Sodium Nitroprusside Reagent Test for Methamphetamine/MDMA
- 1- NARK200VCC Compact Carrying Case with Foam Insert

No. NARK12PAK	NARK II 12-Pak Notebook Kit
No. NARK24PAK	NARK II 24-Pak Notebook Kit
No. NARK200VC	NARK II Patrol Kit

NARK II Presumptive Drug Test Kits

Popular tests packaged in convenient kits. We offer a variety of combination field tests in rugged plastic cases from 16 to 130 individual tests. Custom kits are available.

NARK200CD Contents:

- 10- NARK2001-1 Marquis Reagent (Opium Alkaloids)
- 10- NARK2004-1 Ehrlich's Reagent (LSD)
- 10- NARK2007-1 Scott Reagent Modified (Cocaine HCl, Crack)
- 10- NARK20010-1 Special Opiates Reagent (Fentanyl, Heroin, Oxycodone)
- 10- NARK20015-1 Sodium Nitroprusside Reagent (Methamphetamine, MDMA)
- 1- NARK2006 Acid Neutralizer, 1 oz. (30ml)
- 1- KCP1000 Black Plastic Copolymer Carrying Case; Dimensions: 8.5" x 5.5" x 3.1875" (21.6cm x 14cm x 8.1cm); Weight: 2 lbs (907 g)

NARK200MEG Contents:

- 10- NARK2001-1 Marquis Reagent (Opium Alkaloids)
- 10- NARK2004-1 Ehrlich's Reagent (LSD)
- 10- NARK2007-1 Scott Reagent Modified (Cocaine HCl, Crack)
- 10- NARK20011-1 Mecke's Reagent (All Heroin)
- 10- NARK20015-1 Sodium Nitroprusside Reagent (Methamphetamine, MDMA)
- 1- NARK2006 Acid Neutralizer, 1 oz. (30ml)
- 1- KCP1001 Black Plastic Copolymer Carrying Case; Dimensions: 8.5" x 5.5" x 3.1875" (21.6cm x 14cm x 8.1cm); Weight: 2 lbs (907 g)

NARK200PK Contents:

- 2- NARK2001-1 Marquis Reagent (Opium Alkaloids)
- 2- NARK2002-1 Nitric Acid Reagent (Heroin/Morphine)
- 2- NARK2003-1 Dille-Koppanyi Reagent (Barbiturates)
- 2- NARK2004-1 Ehrlich's Reagent (LSD)
- 2- NARK2005-1 Duquenois-Levine Reagent (Marijuana, Hashish, Hash Oil, THC)
- 2- NARK2007-1 Scott Reagent Modified (Cocaine HCl, Crack)
- 2- NARK20011-1 Mecke's Reagent (All Heroin)
- 2- NARK20015-1 Sodium Nitroprusside Reagent (Methamphetamine, MDMA)
- 1- NARK2006 Acid Neutralizer, 1 oz. (30ml)
- 1- Black Plastic Copolymer Carrying Case; Dimensions: 8.5" x 5.5" x 3.1875" (21.6cm x 14cm x 8.1cm); Weight: 1.45 lbs (657 g)

NARK200MK Contents:

- 2- NARK2001 Marquis Reagent, 10 Tests
- 2- NARK2005 Duqenois-Levine Reagent, 10 Tests
- 1- NARK2002 Nitric Acid Reagent, 10 Tests
- 1- NARK2003 Dille-Koppanyi Reagent, 10 Tests
- 1- NARK2004 LSD Reagent, 10 Tests
- 1- NARK2007 Scott Reagent (modifid), 10 Tests
- 1- NARK2008 Methadone Reagent, 10 Tests
- 1- NARK20010 Special Opiates Reagent, 10 Tests
- 1- NARK20011 Mecke's Modified Reagent, 10 Tests
- 1- NARK20012 Talwin Reagent, 10 Tests
- 1- NARK20015 Sodium Nitroprusside Reagent, 10 Tests
- 1- NARK2006 Acid Neutralizer, 1 oz. (30ml)
- 1- Black Plastic Carrying Case

No. NARK200CD	NARK II Customs/Douanes Kit
No. NARK200MEG	NARK II Meg Kit
No. NARK200PK	NARK II Port-A-Kit
No. NARK200MK	NARK II Master Kit



Mollies Drug Reagent Kit

The Mollies Kit contains 50 individual tests to help you presumptively identify substances commonly known as Mollies on the street.

NARK200MLE Contents:

- 2- NARK20032 Mollies Reagent, 10/box
- 1- NARK20029 2C Reagent, 10/box
- 1- NARK20015 Meth/MDMA (Ecstasy) Reagent, 10/box
- 1- NARK20024 MDPV (Synthetic Cathinones) Reagent, 10/box
- 1- KCP1000 Plastic Copolymer Carrying Case
- 1- Mollies Sequential Testing Chart

No. NARK200MLE	NARK II Mollies Drug Reagent Kit
----------------	----------------------------------

Master Synthetic Drug Reagent Kit

With this kit, you'll have all our reagent kits to conduct a presumptive drug identification for synthetic substances in one easy and convenient package. There's an acid neutralizer included so presumptive tests can be neutralized and you can dispose of the test according to your department protocols.

NARK200SYN Contents:

- 3- NARK20023 Synthetic Cannabinoids Reagent, 5/box
- 2- NARK20032 Mollies Reagent, 10/box
- 2- NARK20015 Methamphetamine/MDMA (Ecstasy), 10/box
- 1- NARK20010 Special Opiates Reagent, 10/box
- 1- NARK20024 MDPV (Synthetic Cathinones) Reagent, 10/box
- 1- NARK20025 Mephedrone (Synthetic Cathinones) Reagent, 10/box
- 1- NARK20026 a-PVP (Synthetic Cathinones) Reagent, 10/box
- 1- NARK20029 2C Reagent, 10/box
- 1- NARK20030 Psilocybin/Psilocin Reagent, 10/box
- 1- NARK2006 Acid Neutralizer Liquid 1 oz. (30ml)
- 1- KCP1002 Black Plastic Copolymer Carrying Case; Dimensions: 8.5"x5.5"x3.1875" (21.6cm x 14cm x 8.1cm); Weight 2 lbs (907g)

No. NARK200SYN	NARK II Master Synthetic Drug Reagent Kit
----------------	---

NARK II Cocaine ID Swipes

Presumptively identify cocaine with only trace amounts left on surfaces. Easy-to-use and results are immediate: the swipe will turn blue if the results are presumptively positive. Sold in 50/box.

No. NS20150	NARK II Cocaine ID Swipes, 50 ea.
-------------	-----------------------------------



NARK II Corrections Kit

Ideal for Corrections Officers, we developed this kit to include the presumptive drug analysis tests for substances commonly abused in prisons. To make it easier for you, we've included a chart to properly identify drugs during tests.

NARK200CR Contents:

- 1- NARK2001 Marquis Reagent (General Screening), 10 Tests
- 2- NARK2005 Duquenois-Levine Reagent (Marijuana, Hash & Hash Oil, THC/ Concentrates), 10 Tests
- 1- NARK2007 Scott Reagent (Cocaine HCl, Base), 10 Tests
- 2- NARK20010 Special Opiates Reagent (Codeine, Heroin, Morphine, Oxycodone, Suboxone), 10 Tests
- 1- NARK20011 Mecke's Modified Reagent (Heroin), 10 Tests
- 1- NARK20015 Sodium Nitroprusside Reagent (Meth, MDMA [Ecstasy]), 10 Tests
- 4- NARK20023 Synthetic Cannabinoids Reagent (Indole formulas only), 5 Tests
- 1- NARK20029 2C Synthetics Reagent (2C Substances including Smiles), 10 Tests
- 1- KCP1002 Copolymer Carry Case
- 1- NARK Corrections Sequential Chart

No. NARK200CR	NARK II Corrections Kit
---------------	-------------------------

NARK II International Kit

This kit was originally designed to appeal to international agencies and comes with a total of 130 tests plus acid nuetralizer.

NARK200IK Contents:

- 10- NARK2001-1 Marquis Reagent (Opium Alkaloids)
- 10- NARK2002-1 Nitric Acid Reagent (Heroin/Morphine)
- 10- NARK2003-1 Dille-Koppanyi Reagent (Barbiturates)
- 10- NARK2004-1 Ehrlich's Reagent (LSD)
- 10- NARK2005-1 Duquenois-Levine Reagent (Marijuana, Hashish, Hash Oil, THC)
- 10- NARK2007-1 Scott Reagent Modified (Cocaine HCl, Crack)
- 10- NARK2008-1 Methadone Reagent
- 10- NARK20010-1 Special Opiates Reagent (Fentanyl, Heroin, Oxycodone)
- 10- NARK20011-1 Mecke's Reagent (All Heroin)
- 10- NARK20012 -1 Talwin (Pentazocine) Reagent
- 10- NARK20013-1 Ephedrine Reagent (Ephedrine, Psuedoephedrine)
- 10- NARK20014-1 Valium, Rohypnol, Ketamine Reagent
- 10- NARK20015-1 Sodium Nitro Prusside Reagent (Methamphetamine, MDMA)
- 1- NARK2006 Acid Neutralizer, 1 oz. (30ml)
- 1- KCP1002 Black Plastic Copolymer Carrying Case; Dimensions: 12.688" x 8.938" x 3.625" (32.22cm x 22.7cm x 9.2cm); Weight: 5.2 lbs (2358 g)

No. NARK200IK	NARK II International Kit
---------------	---------------------------



**Pricing Information:**  
For current pricing visit  
www.sirchie.com, or  
call 800.356.7311

**Often Bought Together:**  
Indicates items often  
purchased in conjunction  
with other items.

NAK Tube Drug Test Reagents

NAK® Reagents are sealed with thin glass ampoules packaged in disposable tubettes—a unique, low cost, easy-to-handle format. Presumptive results are based on a color reaction of the reagents with suspect substances.

**These presumptive results should always be confirmed by an analytical laboratory.**

**Often Bought Together:**  
(SF0081L) Black Nitrile  
Gloves, Large



No.	Description	Presumptively tests for
NAR10001	Mayer's Reagent, 10/box	Narcotic Alkaloids
NAR10002	Marquis Reagent, 10/box	Opium Alkaloids
NAR10003	Nitric Acid Reagent, 10box	Differentiate Heroin/Morphine
NAR10004	Cobalt-Thiocyanate Reagent, 10/box	Cocaine/Crack
NAR10005	Dille-Koppanyi Reagent, 10/box	Barbiturates
NAR10006	Mandelin Reagent, 10/box	Amphetamines
NAR10007	Ehrlich's Reagent, 10/box	LSD
NAR10008	Duquenois Reagent, 10/box	Marijuana, Hash, Hash Oil, THC Concentrates
NAR10009	KN Reagent, 10/box	Marijuana (green plant & seeds)
NAR10013	Modified Cobalt-Thiocyanate Reagent, 10/box	Cocaine Freebase, Crack
NAR10014	Valium Reagent, 10/box	Valium, Rohypnol, Methcathinone, Ketamine
NAR10015	Meth/MDMA Reagent, 10/box	Methamphetamine, MDMA (Ecstasy)
NAR10016	Mecke's Modified Reagent, 10/box	Heroin (white, brown, black tar)
NAR10018	Frohdess Reagent, 10/box	Talwin® (Pentazocine), DMT (Dimethyltryptamine)
NAR10019	Ephedrine Reagent, 10/box	Ephedrine, Pseudoephedrine
NAR10022	Special Opiates Reagent, 10/box	Codeine, Heroin, Oxycodone
NAR10024	MDPV Reagent, 10/box	MDPV (Synthetic Cathinones)
NAR10025	Mephedrone Reagent, 10/box	Mephedrone (Synthetic Cathinones)
NAR10028	GHB Reagent, 5/box	GHB (Gamma Hydroxybutyrate)
NAR10029	PCP/Methaqualone Reagent, 10/box	PCP, Methaqualone
NAR10030	Psilocybin/Psilocin Reagent, 10/box	Psilocybin, Psilocin (mushrooms)
NAR10010	Polyester Bags, 10/box	
NAR10011	Acid Neutralizer, 1 oz. Bottle	

NAK Analysis Reagent Kit

We included reagents for the most frequently abused drugs in this field test kit so you always have them on hand. Each individual reagent is sealed in a thin, glass ampoule which eliminates cross-contamination with other reagents.

- NA100 Contents:**
- 1- NA10001 Mayer's Reagent for Narcotics Alkaloids, Set of 10 Tests

1- NA10002 Marquis Reagent for Opium Alkaloids, Amphetamines/Meth, 10 Tests

1- NA10003 Nitric Acid Reagent to differentiate between Heroin and Morphine, 10 Tests

1- NA10004 Cobalt Thiocyanate Reagent for Cocaine, Crack, 0 Tests

1- NA10005 Dille-Koppanyi Reagent for Barbiturates, 10 Tests

1- NA10006 Mandelin Reagent for Amphetamines, 0 Tests

1- NA10007 Ehrlich's Reagent for LSD, 10 Tests

1- NA10008 Duquenois Reagent for Marijuana, Hashish and THC, 10 Tests

1- NA10010 Poly Bags, 10

1- NA10011 Acid Neutralizer

1- NA10012 Black, Plastic, Copolymer Carrying Case, 8.5" x 5.5" x 3.1875"

No. NAR100

NAK Analysis Reagent Kit

Thin Layer Chromatography of Marijuana Kit

This is a rapid, simplified, specific test for marijuana using a laboratory technique that any officer can perform. Includes six individual tests.

- MZ020 Contents:**
- 6- MZ020A Thin Layer Chromatography Test

1- MZ010B Extract Reagent #1, 1 oz.

1- KCP228A Glass Jar with Closure, 4 oz., Labeled Mobile Phase

1- KCP228A Glass Jar with Closure, 4 Oz, Labeled Color Developing

1- 288DP Disposable Pipettes with Capillary and 3ml Bulb

1- EIL011 Evidence Identification Labels, "Crime Evidence"

1- KCP321 Photo: Chromatography of Marijuana Comparison

1- KCP1391 Plastic Tweezers, 4.75"

1- KCP165 Glass Stirring Rod, 5"

1- KCP154 Pipette with Suction Bulb

1- MZ020C Texturized, Molded-Plastic Carrying Case, 14.25" x 8.25" x 7"

No. MZ020

Thin Layer Chromatography of Marijuana Kit

No. MZ020A

Replacement Test, 1 ea.

Fen-Her™ Fentanyl/Heroin Test Kit

**Detect and differentiate trace or visible amounts of fentanyl and heroin**

Sirchie is proud to offer the Fen-Her™ Kit for law enforcement. The Fen-Her kit is a new, safer and highly accurate kit for the presumptive identification of fentanyl, fentanyl derivatives, and heroin.

Fen-Her™ requires a very small amount of sample. The user collects a sample of suspected particle(s), residues from unopened bags or other surfaces. The detection limit for fentanyl is approximately 135 nanograms (135 billionths of a gram); the detection limit for heroin is approximately 175 nanograms. The entire sampling and analyzing procedure takes just a few minutes.

- Safer because the sample pen needs just a trace amount of material - samples can be taken from residues on paraphernalia such as spoons, pipes, and edges of plastic bags, etc.
- Accurate because color appears in different locations on a test card eliminating confusion between heroin, fentanyl & derivatives
- Separates common cutting agents and adulterants such as diphenhydramine (Benadryl™)



**COLORS** - fentanyl is an orange color; heroin is a purple color. If an orange color appears outside of the fentanyl zone, it isn't fentanyl. If a purple spot appears outside the heroin zone, it isn't heroin.

**ZONES** - fentanyl and heroin appear in zones. The actual position heroin and fentanyl can vary vertically. The ZONES are marked on the inside of WRAPPER to be used as a visual guide.

No. NARK200FENHER

Fen-Her™ Fentanyl/Heroin Test Kit, 5 tests





**Pricing Information:**  
For current pricing visit  
www.sirchie.com, or  
call 800.356.7311

**Often Bought Together:**  
Indicates items often  
purchased in conjunction  
with other items.

NARKsafe™ NDA  
Narcotic Disposal  
Agent

NARKsafe NDA is a drug disposal solution that enables the onsite disposal of drug evidence (20 grams or less) – quickly, easily and safely, rendering the drugs unusable.

The only polymer-based drug disposal solution on the market, NARKsafe NDA enables onsite disposal in an effective, environmentally-safe method that can be tossed into the trash.

- Dispose of drug evidence for post-adjudicated cases that have been dismissed, pled out or where a destructive order has been issued
- Eliminates the cost of holding onto small amounts of adjudicated drug evidence in storage until a quantity large enough to meet minimum destruction requirements at (incinerator) disposal facilities is accumulated
- Reduces the risk for abuse, theft and misuse associated with keeping drugs in evidence storage vaults
- Independent third-party testing includes tests on retrievability, extractability, toxicity characteristics leaching procedure and hazardous waste screening have been conducted to ensure that NARKsafe NDA is effective, safe, environmentally friendly and with no risk of retrievability.
- For use with 20 grams or less of powder, pills, tablets or capsules.

No. NDAK100	NARKsafe NDA Kit, 1 each
-------------	--------------------------

Chemical Splash Goggles

Protect your eyes from chemicals in the lab. Full coverage goggles fit over regular eye glasses.



No. CSG1001	Chemical Splash Goggles
-------------	-------------------------

Black, Powder-Free  
Nitrile Gloves

These thin, comfortable, ambidextrous nitrile gloves provide extra puncture and abrasion resistance. Always use nitrile gloves while working with illicit substances and NARK presumptive drug kits.



No. SF0081S	Black, Powder-Free, Sm, 100 ea.
No. SF0081M	Black, Powder-Free, Md, 100 ea.
No. SF0081L	Black, Powder-Free, Lg, 100 ea.
No. SF0081XL	Black, Powder-Free, XL, 100 ea.

Blue, Powder-Free  
Nitrile Gloves

These blue, powder-free, nitrile gloves provide extra puncture and abrasion resistance over latex. Sold in boxes of 100 ambidextrous gloves.



No. SF00851	Blue, Powder-Free, Sm, 100 ea.
No. SF00852	Blue, Powder-Free, Md, 100 ea.
No. SF00853	Blue, Powder-Free, Lg, 100 ea.
No. SF00854	Blue, Powder-Free, XL, 100 ea.

Long Cuff Nitrile Gloves

These long cuff, powder-free, durable, 8 mil nitrile gloves provide extra thickness for added protection. Perfect for tough jobs where extra durability or safety is needed. Meets NFPA standards for MES workers and provides protection against bloodborne pathogens. Available in five sizes, and come 50 each per box.



No. SF00871	Long Cuff Nitrile Gloves, Sm, 50 ea.
No. SF00872	Long Cuff Nitrile Gloves, Md, 50 ea.
No. SF00873	Long Cuff Nitrile Gloves, Lg, 50 ea.
No. SF00874	Long Cuff Nitrile Gloves, XL, 50 ea.
No. SF00875	Long Cuff Nitrile Gloves, 2XL, 50 ea.

StrongSuit TecArmor Gloves

True cut and puncture resistance packed into a true fit (real anatomic pre-curved palm) for the highest level of protection: CR A5 and PR A5 (ANSI 105). Fully lined with Kevlar®, StrongTip fingers, and padded on the knuckles.



Features:

- **Cut level A6 (ABSI 105):** Designed for the professional who demands the highest level of cut resistance.
- **Puncture level A5 (ANSI 105):** High risk tasks require the highest level of protection. TecArmor Plus keeps you protected from exposure.
- **Fit and Function:** TecArmor gives you high levels of CR and PR protection without sacrificing fit and feel. This glove is pre-curved and styled to give you the best fit and protection available.

No. SG51600S	TecArmor Gloves, small
No. SG51600M	TecArmor Gloves, medium
No. SG51600L	TecArmor Gloves, large
No. SG51600XL	TecArmor Gloves, XL
No. SG51600XXL	TecArmor Gloves, XXL

NARKsafe™ Single-Use  
Powder Substance Packaging,  
Transport, Testing Kit

Any powdery substance encountered on the street should be considered dangerous. Law enforcement needs to proceed with caution while handling, transporting, and testing these substances.

Following DEA and CDC guidelines, we developed this one-time use PPE Kit to be used specifically around powdery substances that could contain fentanyl, carfentanil, or other powerful opioid-based substances.



NARK300PPE Contents:

- 1- Coveralls (choose size)

1- P-100 Respirator Mask

1- Black Nitrile Gloves, pair

1- Blue, Long Cuff nitrile Gloves, pair

1- Boot Covers, pair

1- Fully-sealed Goggles

2- Fentanyl Warning Labels
- 1- Polyester Evidence Integrity bag

1- Plastic zip-top evidence bag

1- Yellow hazardous material bag

1- donning/doffing instruction card for properly using personal protection gear

No. NARK300PPSM	NARKsafe™ Single-Use Powder Substance Kit, Md
No. NARK300PPEL	NARKsafe™ Single-Use Powder Substance Kit, Lg
No. NARK300PPEXL	NARKsafe™ Single-Use Powder Substance Kit, XL
No. NARK300PPE2XL	NARKsafe™ Single-Use Powder Substance Kit, 2XL
No. NARK300PPE3XL	NARKsafe™ Single-Use Powder Substance Kit, 3XL

NARKsafe™ Patrol Kit

NARKsafe™ Patrol Kit keeps you safe when you encounter narcotics in the field. The kit includes one respirator mask with 100% filtration of harmful substances, fully sealed goggles and a pair of nitrile gloves. To stay safe, take the same precautions on patrol as they do in the lab, whether conducting presumptive field tests or collecting drug evidence, wear protective gear. ***If fentanyl is suspected, use our NARK300PPE kit for more complete protection.***



NARK300PK Contents:

- 1- P-100 Respirator Mask

1- Fully-sealed Goggles

1- Black Nitrile Gloves, pair

No. NARK300PK	NARKsafe Patrol Kit
---------------	---------------------

Disposable P100  
Respirator Mask

This Respirator offers the highest level of NIOSH rated filtration efficiency in a maintenance-free respirator. It features a 3M™ Cool Flow™ Exhalation Valve, a faceseal, adjustable noseclip, lightweight construction and an all filter design. Sold individually. Meets 42 CFR84 standard.



No. NARK300P100R	Disposable P100 Respirator Mask
------------------	---------------------------------



**Pricing Information:**  
For current pricing visit  
www.sirchie.com, or  
call 800.356.7311



**Often Bought Together:**  
Indicates items often  
purchased in conjunction  
with other items.

Micromax® NS Coveralls

According to the Drug Enforcement Agency (DEA) and Centers for Disease Control (CDC), a full coverall with hood is recommended for handling dangerous illicit substances. Below are our most popular sizes.



**Often Bought Together:**  
Micromax NS  
Shoe & Boot Covers



No. MMH300M	Coveralls, Zipper, Attached Hood, Elastic Wrist & Ankle, Md
No. MMH300L	Coveralls, Zipper, Attached Hood, Elastic Wrist & Ankle, L
No. MMH300XL	Coveralls, Zipper, Attached Hood, Elastic Wrist & Ankle, XL
No. MMH3002XL	Coveralls, Zipper, Attached Hood, Elastic Wrist & Ankle, 2XL

Micromax® NS Boot Covers

Breathable for comfort, these boot covers remain strong whether they're wet or dry.



No. BC17LXL	Boot Covers, Elastic Top, L/XL, pair
No. BC172XL	Boot Covers, Elastic Top, 2XL, pair

Disposable Hazardous Material Bags

Properly dispose of your used personal PPE in our 18" x 30 1/2" hazardous material bags: easily identifiable by their bright yellow color. Always dispose of any hazardous waste according to Federal, state, and local regulations. Sold in packs of 10 bags each.



**Often Bought Together:**  
(NARK300FENL) Fentanyl  
Warning lbel

No. NARK300HZB	Disposable Hazard Bags, 18" x 30 1/2", 10 ea.
----------------	---

NARKsafe™ Packaging Kits

We've put together evidence packaging kits for suspected dangerous substances which include: five (5) Keepsafe Poly Evidence Integrity Bags; five (5) zip-top bags; and five (5) "Danger: Suspect Fentanyl" labels to warn officials handling the bags. Three kit sizes are available: select the large size for packaging over-sized evidence.



No. NARK300PKG01	NARKsafe Packaging Kit 5-pack, Sm
No. NARK300PKG02	NARKsafe Packaging Kit 5-pack, Md
No. NARK300PKG03	NARKsafe Packaging Kit 5-pack, Lg

KeepSafe Poly Evidence Bags

These bags feature a high degree of security: once closed, unauthorized attempts to break the seal will cause "STOP" signs to appear in the orange strips of the tape. Also includes a tear-away tab for evidence tracking.



No. KSB610C	KeepSafe Bags, 6" x 10 1/4", 100 ea.
No. KSB1013C	KeepSafe Bags, 10" x 13", 100 ea.
No. KSB1216C	KeepSafe Bags, 12" x 16 7/8", 100 ea.

Write-Block Zip-Top Evidence Bag

These zip-top bags are made of 0.004" thick polyethylene and feature a write-on block for easy labeling of contents.



No. ZWB02	Zip-Top Write Block Evidence Bag, 9" x 12", 100 ea.
-----------	---

Fentanyl Warning Labels

These attention-getting labels are printed in full color on durable crack-n-peel material, and measure 2.875" x 3.75".



No. NARK300FENL	Fentanyl Warning Labels, 100/pk
-----------------	---------------------------------

Kapak® Clear Heat Seal Polyester Pouches

Eliminate cross-contamination with these clear, strong, heat-sealable pouches; suitable for sharp objects and illegal or abused substances. Available in 2.5 mil and 4.5 mil thicknesses and in a wide variety of sizes.



No. K2B46C	Kapak® Polyester Pouch, 2.5 mil, 4" x 6", 100 ea.
No. K2B46M	Kapak® Polyester Pouch, 2.5 mil, 4" x 6", 2000 ea.
No. K2B68C	Kapak® Polyester Pouch, 2.5 mil, 6 1/2" x 8", 75 ea.
No. K2B68M2	Kapak® Polyester Pouch, 2.5 mil, 6 1/2" x 8", 2000 ea.
No. K2B812L	Kapak® Polyester Pouch, 2.5 mil, 8" x 12", 40 ea.
No. K2B812M2	Kapak® Polyester Pouch, 2.5 mil, 8" x 12", 2000 ea.
No. K2B1012L	Kapak® Polyester Pouch, 2.5 mil, 10" x 12", 30 ea.
No. K2B916L	Kapak® Polyester Pouch, 2.5 mil, 9 1/2" x 16", 20 ea.
No. K2B1216L	Kapak® Polyester Pouch, 2.5 mil, 12" x 16", 15 ea.
No. K4B46C	Kapak® Polyester Pouch, 4.5 mil, 4" x 6", 100 ea.
No. K4B46M	Kapak® Polyester Pouch, 4.5 mil, 4" x 6", 1000 ea.
No. K4B68C	Kapak® Polyester Pouch, 4.5 mil, 6 1/2" x 8", 75 ea.
No. K4B68M	Kapak® Polyester Pouch, 4.5 mil, 6 1/2" x 8", 1000 ea.
No. K4B89L	Kapak® Polyester Pouch, 4.5 mil, 8" x 9 1/2", 50 ea.
No. K4B89M	Kapak® Polyester Pouch, 4.5 mil, 8" x 9 1/2", 1000 ea.
No. K4B812L	Kapak® Polyester Pouch, 4.5 mil, 8" x 12", 40 ea.
No. K4B812M	Kapak® Polyester Pouch, 4.5 mil, 8" x 12", 1000 ea.
No. K4B1012L	Kapak® Polyester Pouch, 4.5 mil, 10" x 12", 30 ea.
No. K4B1012M	Kapak® Polyester Pouch, 4.5 mil, 10" x 12", 1000 ea.
No. K4B916L	Kapak® Polyester Pouch, 4.5 mil, 9 1/2" x 16", 20 ea.
No. K4B9165C	Kapak® Polyester Pouch, 4.5 mil, 9 1/2" x 16", 500 ea.
No. K4B1216L	Kapak® Polyester Pouch, 4.5 mil, 12" x 16", 15 ea.
No. K4B12165C	Kapak® Polyester Pouch, 4.5 mil, 12" x 16", 500 ea.

Heat Sealers for Polyester Pouches and Tubing

Seal the Kapak polyester pouches and tubing with this heat sealer. Available in 10 1/2" and 18" sizes.



Not intended for use with polyethylene plastic bags.

No. BPHS10	Poly Pouch Heat Sealer, 10 1/2" (26.7 cm)
No. BPHS18	Poly Pouch Heat Sealer, 18" (45.7 cm)

Biohazard Evidence

Sharps Container

Ideal for used medical needles or other sharp medical instruments, you should have this container on hand to ensure a safe transport. Available in 1 quart.



No. SHRP32	Sharps Container, 1 qt.
------------	-------------------------

Sharps Collector

These 5.4 quart collectors have a visible full level indicator so they can be safely handled during transport.



No. SHRP54	Sharps Collector, 5.4 qt.
------------	---------------------------

Sharpstube™ Syringe Transport Tubes

Designed to hold a hypodermic syringe for transportation, these tubes are well-marked with bright orange-red international biohazard markings. Sold in pack of 12 each.



No. ECT2	Syringe Transport Tube, 12 ea.
----------	--------------------------------



**Pricing Information:**  
For current pricing visit  
www.sirchie.com, or  
call 800.356.7311

**Often Bought Together:**  
Indicates items often  
purchased in conjunction  
with other items.

## Evidence Collection Tubes

Our clear evidence collection tubes will keep evidence with sharp edges visible and secure during transport.

### Evidence Collection Tubes

Available in two sizes for single syringes or knives, or multiple syringes. An ethafoam stopper will immobilize the objects while they're in transport.



No. ECT1	Evidence Collection Tubes, 1 5/16" x 8", 12 ea.
No. ECT3	Evidence Collection Tubes, 3" x 12", 8 ea.

### Telescopic Evidence Collection Tubes

These tubes have length extensions and a twist-lock mechanism that keeps the tubes in the extended position.



**Often Bought Together:**  
(SM1000) Sirchmark  
Evidence Integrity Tape

No. ECT4	Twist Tube, 1.26" dia. x 4 3/4 - 7 7/8" length, 12 ea.
No. ECT5	Twist Tube, 2.05" dia. x 7 7/8 - 13 3/4" length, 12 ea.
No. ECT6	Twist Tube, 2.56" dia. x 13 3/4 - 24 3/8" length, 12 ea.

## Evidence Integrity Tape

You should use integrity tapes during your narcotics investigation to show a chain of custody and that there hasn't been tampering with evidence. Regular use of these tapes will ensure the integrity of your investigation.

### Sirchmark™ Evidence Tape

This 1 3/8" x 108ft evidence integrity tape features a matte acetate base, a quick-setting acrylic adhesive, and a dual-serrated edge. Tape with a white stripe is ideal for writing information with a ball-point pen. Comes in a handy dispenser box.



No. SM10002	Red, "Evidence", 108ft
No. SM10002NE	Red, "Narcotics Evidence", 108ft
No. SM50002	Red/White Stripe, "Evidence", 108ft
No. SM50002NE	Red/White Stripe, "Narcotics Evidence", 108ft

### Sirchmark™ Evidence Strips

These Evidence Integrity Strips offer the same features as our Evidence Integrity tape, but are cut into 1 3/8" wide x 7" long strips.



No. SM1000SR	Red, "Evidence", 100 ea.
No. SM1000SNE	Red, "Narcotics Evidence", 100 ea.
No. SM1000SRW	Red/White Stripe, "Evidence", 100 ea.
No. SM1000SRWNE	Red/White Stripe, "Narcotics Evidence", 100 ea.

### Basic Compact Scales

These Basic Compact Scales offer all of the necessary features for quick and easy weighing. Lightweight and portable these scales feature battery as well as AC operation, are splashproof and offer 3 weighing units (g, lb, oz).



### Advanced Scales

The Advanced Scale is the best choice for simple operation and economy, featuring a stainless steel pan, three-point overload protection and battery as well as AC power operation.

### Precision Scales

Precision balances have what it takes for basic lab work, field use and various industrial applications. These balances feature 15 weighing units, three-point overload protection, RS-232 and USB interface, internal and external calibration and battery as well as AC operation.

No. SCL200	Basic Compact Scale (200g), 110V
No. SCL200E	Basic Compact Scale (200g), 220V
No. SCL5000	Basic Compact Scale (5000g), 110V
No. SCL5000E	Basic Compact Scale (5000g), 220V
No. SCL2600	Advanced Scale (2600g), 110V
No. SCL2600E	Advanced Scale (2600g), 220V
No. SCL1000	Precision Balance (1000g), 110V
No. SCL1000E	Precision Balance (1000g), 220V

## Liquid Sampling Kits

Many protocols for the collection, transport, and storage of hazardous substances, including those often found at clandestine meth labs, recommend or require placing the substance in a glass jar and then placing the jar inside a Nalgene® jar. The use of leak proof glass and Nalgene® jars helps to enhance the safety of all personnel and preserve the evidence. Each of the Sirchie Liquid Sampling Kits is designed to have a glass jar that fits perfectly into a Nalgene® jar with pipettes for collection of the substance.

### 30ml Liquid Sampling Kit

This kit includes 30ml Nalgene® Bottles (12 each), 4ml Glass Vials (12 each), and Disposable Pipettes (12 each).



No. LSK100	30ml Liquid Sampling Kit
------------	--------------------------

### 60ml Liquid Sampling Kit

This kit includes 60ml Nalgene® Bottles (12 each), 12ml Glass Vials (12 each), and Disposable Pipettes (12 each).



No. LSK200	60ml Liquid Sampling Kit
------------	--------------------------

### 125ml Liquid Sampling Kit

This kit includes 125ml Nalgene® Bottles (12 each), 16ml Glass Vials (12 each), and Disposable Pipettes (12 each).



No. LSK300	125ml Liquid Sampling Kit
------------	---------------------------

### 250ml Liquid Sampling Kit

This kit includes 250ml Nalgene® Bottles (12 each), 30ml Glass Vials (12 each), and Disposable Pipettes (12 each).



No. LSK400	250ml Liquid Sampling Kit
------------	---------------------------

### Disposable Liquid Sampling Pipettes

Use these disposable pipettes for sampling unknown liquids. Store collected liquids in glass vials placed inside Nalgene bottles.



No. PIP25	Disposable Liquid Sampling Pipettes, 25 ea.
No. PIP500	Disposable Liquid Sampling Pipettes, 500 each

### Translucent Nalgene® Bottles

These bottles can be used to perform general purpose laboratory applications. These bottles are leakproof, light-weight, and a more flexible version of Nalgene® lab quality bottles. Translucent, natural colored polypropylene screw cap included with each bottle.



HDPE containers can be used for long-term storage but direct UV exposure should be avoided

#### Features:

- Economical and recyclable
- Wide mouth makes these bottles easy to fill with solids or liquids
- Linerless caps ensure leakproof performance without a liner that can wrinkle, cause leaks and contaminate your solutions

No. NAJAR0001	Translucent Nalgene® Bottles, 30ml, 12 ea.
No. NAJAR0002	Translucent Nalgene® Bottles, 60ml, 12 ea.
No. NAJAR0004	Translucent Nalgene® Bottles, 125ml, 12 ea.
No. NAJAR0008	Translucent Nalgene® Bottles, 250ml, 12 ea.
No. NAJAR0016	Translucent Nalgene® Bottles, 500ml, 12 ea.
No. NAJAR0032	Translucent Nalgene® Bottles, 1 liter, 6 ea.

### Amber Nalgene® Lab Quality Bottles

Store, ship and package light-sensitive liquid using these narrow-mouth Amber HDPE lab quality Nalgene® bottles. Featuring excellent chemical resistance to most acids, bases, and alcohols, these Nalgene® bottles are high density polyethylene and are guaranteed leakproof. Includes amber polypropylene screw caps for each bottle.



#### Features:

- Opaque amber reduces UV light transmissions to protect light-sensitive liquids
- Meets U.S. Pharmacopoeia 23 requirements for light transmission
- Versatile bottles where exposure to light needs to be restricted
- Can be used in freezer to -100°C (-148°F) for long term storage
- Linerless caps provide the ultimate in leakproof protection without the use of liner that can wrinkle, cause leaks or contaminate your reagents

No. NABOT9025	Amber Nalgene® Bottles, 8ml, 12 ea.
No. NABOT9050	Amber Nalgene® Bottles, 15ml, 150 ea.
No. NABOT0001	Amber Nalgene® Bottles, 30ml, 12 ea.
No. NABOT0002	Amber Nalgene® Bottles, 60ml, 12 ea.
No. NABOT0008	Amber Nalgene® Bottles, 250ml, 12 ea.
No. NABOT0016	Amber Nalgene® Bottles, 500ml, 12 ea.
No. NABOT0032	Amber Nalgene® Bottles, 1 liter, 6 ea.

**Pricing Information:**  
For current pricing visit  
[www.sirchie.com](http://www.sirchie.com), or  
call 800.356.7311



**Why NARKsafe™ is the only enclosure to use when testing illicit substances:**

- Horizontal laminar airflow pattern
- Dual HEPA filtration
- HEPA filters with 99.997% efficiency rating for collecting 0.3-micron particles
- Electrostatically charged pre-filter to collect larger particulate and minimize clogging of the HEPA filter
- ASHRAE 110 testing adjustable face velocity
- Ability to replace contaminated filters under negative pressure
- Blue flooring to see substances clearly when testing




When testing illicit substances and changing filters, the operator should always follow Drug Enforcement Administration (DEA) and the Centers for Disease Control (CDC) guidelines for handling fentanyl, carfentanil and acetyl-fentanyl.

**NARKsafe™ Illicit Substance Testing Workstation**

Fentanyl-related deaths are increasing as the opioid epidemic increases worldwide. Fentanyl is deadly if swallowed, inhaled, or absorbed through the skin, so being protected when handling suspected substances is paramount. NARKsafe™ fume hood enclosures are ideal for protecting professionals during the forensic analysis of fentanyl and other dangerous substances. Only this system offers a turbulence-free, horizontal laminar airflow with dual HEPA filtration.



NARKsafe™ enclosure provides a negative pressure environment to pull any particulate away from the operator and eliminate end-user exposure. It also allows for the addition of a bonded carbon filter, so analysis of other chemical solutions or marijuana exhibits can be safely handled in the same workspace. Available in three widths: 36", 48", 72".



No. ACDD760	NARKsafe™ 36" Powder Safety Workstation, 110V AC w/prefilter, carbon filter & HEPA filters
No. ACDD760220	NARKsafe™ 36" Powder Safety Workstation, 220V AC w/prefilter, carbon filter & HEPA filters
No. ACDD770	NARKsafe™ 48" Powder Safety Workstation, 110V AC w/prefilter, carbon filter & HEPA filters
No. ACFPRE6	Replacement Prefilter for NARKsafe, 12 pack
No. ACEHEPA18	Replacement 1st HEPA Filter for NARKsafe, 1 ea.
No. ACFHEPA	Replacement 2nd HEPA Filter for NARKsafe, 1 ea.
No. ACF100	Replacement Carbon Filter for NARKsafe, 1 ea.



Featured above: Sterile round foam-tipped applicators (FTS50W).



From evidence collection to gun shot residue collection to DNA collection to sexual assault evidence collection, we can customize your kit and needs. We work closely with states and jurisdictions to ensure their kits have all the necessary components to allow for a thorough and easy examination.

(800) 356-7311 | [sales@sirchie.com](mailto:sales@sirchie.com) | [www.sirchie.com](http://www.sirchie.com)