



# Narcotics Investigation

As the opioid epidemic grows and becomes more dangerous to law enforcement and first responders, our goal is to make you as safe as possible. We've reviewed the Drug Enforcement Administration (DEA) and the Centers for Disease Control (CDC) guidelines for handling fentanyl, carfentanil and acetyl-fentanyl.

Building from guidelines offered by the Drug Enforcement Administration (DEA), the Centers for Disease Control (CDC), and the manufacturers of pharmaceutical fentanyl, we offer solutions under our NARKsafe™ product line.

NARKsafe™ solutions are designed to help keep you and your team safe while handling, packaging, transporting, testing, and storing potential fentanyl, carfentanil, acetyl-fentanyl, and other dangerous substances.

**For more information about these guidelines, please see:**

- Fentanyl: A Briefing Guide for First Responders, DEA [https://www.dea.gov/druginfo/Fentanyl\\_BriefingGuideforFirstResponders\\_June2017.pdf](https://www.dea.gov/druginfo/Fentanyl_BriefingGuideforFirstResponders_June2017.pdf)
- Fentanyl: Preventing Occupational Exposure to Emergency Responders, CDC <https://www.cdc.gov/niosh/topics/fentanyl/default.html>

**NARKsafe™**



## Fentanyl Facts:



- Fentanyl is a synthetic opioid developed as pain management treatment
- 50 to 100 times more potent than morphine and 25 to 40 times more potent than heroin
- Appearance: usually white powder similar to methamphetamine and cocaine
- Fentanyl is a Schedule II substance (Heroin is Schedule I)
- Attractive to drug abusers due to its intense high
- Street names: China Girl, China White, Goodfella, Apache, TNT, Tango & Cash, Dance Fever, Jackpot
- Trade Names: Actiq®, Fentora™, Duragesic®, Sublimaze®
- Clandestinely produced and either cut into heroin or mixed with other adulterants and sold as heroin
- Signs of opioid overdose: pinpoint pupils, respiratory depression and unconsciousness
- Antidote: Narcan (Naloxone) availability: <https://www.bjatraining.org/tools/naloxone/Naloxone-Background>

## Carfentanil Facts:

- 100 times more powerful than fentanyl
- Medically used to sedate large animals such as buffalo and elephants
- Two milligrams of carfentanil can sedate a 2,000 pound African elephant

## NARKsafe™ Single-Use Powder Substance Packaging, Transport, Testing Kit

Any powdery substance encountered on the street should be considered dangerous. Law enforcement needs to proceed with caution while handling, transporting, and testing these substances.

Following DEA and CDC guidelines, we developed this one-time use PPE Kit to be used specifically around powdery substances that could contain fentanyl, carfentanyl, or other powerful opioid-based substances.

Kit contains: coveralls; P-100 respirator; two pair of nitrile gloves; goggles; evidence bag; a biohazard bag; and donning/doffing instruction card.



No. NARK300PPEM	NARKsafe™ Single-Use Powder Substance Kit, Md
No. NARK300PPEL	NARKsafe™ Single-Use Powder Substance Kit, Lg
No. NARK300PPEXL	NARKsafe™ Single-Use Powder Substance Kit, XL
No. NARK300PPE2XL	NARKsafe™ Single-Use Powder Substance Kit, 2XL
No. NARK300PPE3XL	NARKsafe™ Single-Use Powder Substance Kit, 3XL

## Disposable P100 Respirator Mask

This Respirator offers the highest level of NIOSH rated filtration efficiency in a maintenance-free respirator. It features a 3M™ Cool Flow™ Exhalation Valve, a face seal, adjustable noseclip, lightweight construction and an all filter design. Sold individually. Meets 42 CFR84 standard.



No. NARK300P100R	Disposable P100 Respirator Mask
------------------	---------------------------------

## Chemical Splash Goggles

Protect your eyes from chemicals in the lab. Full coverage goggles fit over regular eye glasses.



No. CSG1001

Chemical Splash Goggles

## Black Powder-Free Nitrile Gloves

These comfortable, ambidextrous nitrile gloves provide extra puncture and abrasion resistance. Always use nitrile gloves while working with illicit substances and NARK presumptive drug kits.



No. SF0081S	Black, Powder-Free, Sm, 100 ea.
No. SF0081M	Black, Powder-Free, Md, 100 ea.
No. SF0081L	Black, Powder-Free, Lg, 100 ea.
No. SF0081XL	Black, Powder-Free, XL, 100 ea.

## Blue, Powder-Free Nitrile Gloves

These blue, powder-free, nitrile gloves provide extra puncture and abrasion resistance over latex. Sold in boxes of 100 ambidextrous gloves.



No. SF00851	Blue, Powder-Free, Sm, 100 ea.
No. SF00852	Blue, Powder-Free, Md, 100 ea.
No. SF00853	Blue, Powder-Free, Lg, 100 ea.
No. SF00854	Blue, Powder-Free, XL, 100 ea.

## Long Cuff Nitrile Gloves

These long cuff, powder-free, durable, 8 mil nitrile gloves provide extra thickness for added protection. Perfect for tough jobs where extra durability or safety is needed. Meets NFPA standards for MES workers and provides protection against bloodborne pathogens.



No. SF00871	Long Cuff Nitrile Gloves, Sm, 50 ea.
No. SF00872	Long Cuff Nitrile Gloves, Md, 50 ea.
No. SF00873	Long Cuff Nitrile Gloves, Lg, 50 ea.
No. SF00874	Long Cuff Nitrile Gloves, XL, 50 ea.
No. SF00875	Long Cuff Nitrile Gloves, 2XL, 50 ea.



## Micromax® NS Coveralls

According to the Drug Enforcement Agency (DEA) and Centers for Disease Control (CDC), a full coverall with hood is recommended for handling dangerous illicit substances. Below are our most popular sizes.



Often bought with:  
Micromax NS  
Shoe & Boot Covers



No. MMH300M	Coveralls, Zipper, Attached Hood, Elastic Wrist & Ankle, Md
No. MMH300L	Coveralls, Zipper, Attached Hood, Elastic Wrist & Ankle, L
No. MMH300XL	Coveralls, Zipper, Attached Hood, Elastic Wrist & Ankle, XL
No. MMH3002XL	Coveralls, Zipper, Attached Hood, Elastic Wrist & Ankle, 2XL

## Micromax® NS Boot Covers

Breathable for comfort, they remain strong whether they're wet or dry.



No. BC17SM	Boot Covers, Elastic Top, Sm/Md, pair
No. BC17LXL	Boot Covers, Elastic Top, L/XL, pair
No. BC172XL	Boot Covers, Elastic Top, 2XL, pair

## Disposable Hazardous Material Bags

Properly dispose of your used personal PPE in our 18" x 30 1/2" hazardous material bags: easily identifiable by its bright yellow color. Always dispose of any hazardous waste according to Federal, state, and local regulations. Sold in packs of 10 bags each.



No. NARK300H2B	Disposable Hazard Bags, 18" x 30 1/2", 10 ea.
----------------	---

## NARKsafe™ Packaging Kits


We've put together evidence packaging kits for suspected dangerous substances which include: five (5) Keepsafe Poly Evidence Integrity Bags; five (5) zip-top bags; and five (10) "Danger: Suspect Fentanyl" labels to warn officials handling the bags. Three kit sizes are available: select the large size for packaging over-sized evidence.



No. NARK300PKG01	NARKsafe Packaging Kit 5-pack, Sm
No. NARK300PKG02	NARKsafe Packaging Kit 5-pack, Md
No. NARK300PKG03	NARKsafe Packaging Kit 5-pack, Lg

## KeepSafe Poly Evidence Bags

These bags feature a high degree of security: once closed, unauthorized attempts to break the seal will cause "STOP" signs to appear in the orange strips of the tape. Also includes a tear-away tab for evidence tracking.



No. KSB610C	KeepSafe Bags, 6" x 10 1/4", 100 ea.
No. KSB1013C	KeepSafe Bags, 10" x 13", 100 ea.
No. KSB1216C	KeepSafe Bags, 12 x 16 7/8", 100 ea.

## Write-Block Zip-Top Evidence Bag

These zip-top bags are made of 0.004" thick polyethylene and feature a write-on block for easy labeling of contents.



No. ZWB02	Zip-Top Write Block Evidence Bag, 9" x 12", 100 ea.
-----------	---

## Fentanyl Warning Labels

These attention-getting labels are printed in full color on durable crack-n-peel material, and measure 2.875" x 3.75".



No. NARK300FENL	Fentanyl Warning Labels, 100/pk
-----------------	---------------------------------

## Kapak® Clear Heat Seal Polyester Pouches

Eliminate cross-contamination with these clear, strong, heat-sealable pouches; suitable for sharp objects and illegal or abused substances. Available in 2.5 mil and 4.5 mil and in sizes 4" x 6", 6 1/2" x 8", 8" x 9 1/2" (only 4.5mil), 8" x 12", 10" x 12", 9 1/2" x 16", 12" x 16".



**Often Bought with:**  
(BPHS10, BPHS18) Polyester Pouch  
Heat Sealers

No. K2B46C	Kapak® Polyester Pouch, 2.5 mil, 4" x 6", 100 ea.
No. K2B46M	Kapak® Polyester Pouch, 2.5 mil, 4" x 6", 2000 ea.
No. K2B68C	Kapak® Polyester Pouch, 2.5 mil, 6 1/2" x 8", 75 ea.
No. K2B68M2	Kapak® Polyester Pouch, 2.5 mil, 6 1/2" x 8", 2000 ea.
No. K2B812L	Kapak® Polyester Pouch, 2.5 mil, 8" x 12", 40 ea.
No. K2B812M2	Kapak® Polyester Pouch, 2.5 mil, 8" x 12", 2000 ea.
No. K2B1012L	Kapak® Polyester Pouch, 2.5 mil, 10" x 12", 30 ea.
No. K2B1012M2	Kapak® Polyester Pouch, 2.5 mil, 10" x 12", 2000 ea.
No. K2B916L	Kapak® Polyester Pouch, 2.5 mil, 9 1/2" x 16", 20 ea.
No. K2B1216L	Kapak® Polyester Pouch, 2.5 mil, 12" x 16", 15 ea.
No. K4B46C	Kapak® Polyester Pouch, 4.5 mil, 4" x 6", 100 ea.
No. K4B46M	Kapak® Polyester Pouch, 4.5 mil, 4" x 6", 1000 ea.
No. K4B68C	Kapak® Polyester Pouch, 4.5 mil, 6 1/2" x 8", 75 ea.
No. K4B68M	Kapak® Polyester Pouch, 4.5 mil, 6 1/2" x 8", 1000 ea.
No. K4B89L	Kapak® Polyester Pouch, 4.5 mil, 8" x 9 1/2", 50 ea.
No. K4B89M	Kapak® Polyester Pouch, 4.5 mil, 8" x 9 1/2", 1000 ea.
No. K4B812L	Kapak® Polyester Pouch, 4.5 mil, 8" x 12", 40 ea.
No. K4B812M	Kapak® Polyester Pouch, 4.5 mil, 8" x 12", 1000 ea.
No. K4B1012L	Kapak® Polyester Pouch, 4.5 mil, 10" x 12", 30 ea.
No. K4B1012M	Kapak® Polyester Pouch, 4.5 mil, 10" x 12", 1000 ea.
No. K4B916L	Kapak® Polyester Pouch, 4.5 mil, 9 1/2" x 16", 20 ea.
No. K4B9165C	Kapak® Polyester Pouch, 4.5 mil, 9 1/2" x 16", 500 ea.
No. K4B1216L	Kapak® Polyester Pouch, 4.5 mil, 12" x 16", 15 ea.
No. K4B12165C	Kapak® Polyester Pouch, 4.5 mil, 12" x 16", 500 ea.

## Biohazard Evidence

### Sharps Container

Ideal for used medical needles or other sharp medical instruments, you should have this container on hand to ensure a safe transport. Available in 1 quart.



No. SHRP32

Sharps Container, 1 qt.

### Sharps Collector

These 5.4 quart collectors have a visible full level indicator so they can be safely handled during transport.



No. SHRP54

Sharps Collector, 5.4 qt.

### Sharpstube™ Syringe Transport Tubes

Designed to hold a hypodermic syringe for transportation, these tubes are well-marked with bright orange-red international biohazard markings. Sold in pack of 12 each.



No. ECT2

Syringe Transport Tube, 12 ea.



## Evidence Collection Tubes

Our clear evidence collection tubes will keep evidence with sharp edges visible and secure during transport.

### Evidence Collection Tubes

Available in two sizes for single syringes or knives, or multiple syringes. An ethafoam stopper will immobilize the objects while they're in transport.



**No. ECT1** Evidence Collection Tubes, 1 5/16" x 8", 12 ea.

**No. ECT3** Evidence Collection Tubes, 3" x 12", 8 ea.

### Telescopic Evidence Collection Tubes

These tubes have length extensions and a twist-lock mechanism that keeps the tubes in the extended position.



**Often Bought with:**  
(SM1000) Sirchie  
Evidence Integrity Tape

**No. ECT4** Twist Tube, 1.26" dia. x 4 3/4 - 7 7/8" length, 12 ea.

**No. ECT5** Twist Tube, 2.05" dia. x 7 7/8 - 13 3/4" length, 12 ea.

**No. ECT6** Twist Tube, 2.56" dia. x 13 3/4 - 24 3/8" length, 12 ea.

## Evidence Integrity Tape

You should use integrity tapes during your narcotics investigation to show a chain of custody and that there hasn't been tampering with evidence. Regular use of these tapes will ensure the integrity of your investigation.

### Sirchie™ Evidence Tape

This 1 3/8" x 108ft evidence integrity tape features a matte acetate base, a quick-setting acrylic adhesive, and a dual-serrated edge. Tape with a white stripe is ideal for writing information with a ball-point pen. Comes in a handy dispenser box.



**No. SM10002** Red, "Evidence", 108ft

**No. SM10002NE** Red, "Narcotics Evidence", 108ft

**No. SM50002** Red/White Stripe, "Evidence", 108ft

**No. SM50002NE** Red/White Stripe, "Narcotics Evidence", 108ft

### Sirchie™ Evidence Strips

These Evidence Integrity Strips offer the same features as our Evidence Integrity tape, but are cut into 1 3/8" wide x 7" long strips.



**No. SM1000SR** Red, "Evidence", 100 ea.

**No. SM1000SNE** Red, "Narcotics Evidence", 100 ea.

**No. SM1000SRW** Red/White Stripe, "Evidence", 100 ea.

**No. SM1000SRWNE** Red/White Stripe, "Narcotics Evidence", 100 ea.

### TurtleSkin® Duty Gloves

TurtleSkin Duty Gloves give you excellent safety and protection from accidental needle stick, punctures, and cuts to the palm side of the hand, fingers, and finger tips.



**No. TSK006DS** TurtleSkin® Duty Gloves, Small, pair

**No. TSK006DM** TurtleSkin® Duty Gloves, Medium, pair

**No. TSK006DL** TurtleSkin® Duty Gloves, Large, pair

**No. TSK006DXL** TurtleSkin® Duty Gloves, XL, pair

**No. TSK006D2XL** TurtleSkin® Duty Gloves, 2XL, pair

### TurtleSkin® Bravo Gloves

Made with a rugged, mesh knit, these gloves have good breathability and cut and hypodermic needle puncture protection. Very lightweight, this glove is perfect for law enforcement officials in warm weather.



**No. TSK009BS** TurtleSkin® Bravo Gloves, Small, pair

**No. TSK009BM** TurtleSkin® Bravo Gloves, Medium, pair

**No. TSK009BL** TurtleSkin® Bravo Gloves, Large, pair

**No. TSK009BXL** TurtleSkin® Bravo Gloves, XL, pair

**No. TSK009B2XL** TurtleSkin® Bravo Gloves, 2XL, pair

### TurtleSkin® Alpha Gloves

The best selling TurtleSkin Alpha Glove gives you excellent cut and hypodermic needle protection. The palm gives you both cut and needlestick puncture resistance. TurtleSkin Alpha Gloves offer a high level of dexterity and tactile sensitivity.



**No. TSK008AS** TurtleSkin® Alpha Gloves, Small, pair

**No. TSK008AM** TurtleSkin® Alpha Gloves, Medium, pair

**No. TSK008AL** TurtleSkin® Alpha Gloves, Large, pair

**No. TSK008AXL** TurtleSkin® Alpha Gloves, XL, pair

**No. TSK008A2XL** TurtleSkin® Alpha Gloves, 2XL, pair





**Often Bought Together:**  
Indicates items often  
purchased in conjunction with  
other items.

## Liquid Sampling Kits

Many protocols for the collection, transport, and storage of hazardous substances, including those often found at clandestine meth labs, recommend or require placing the substance in a glass jar and then placing the jar inside a Nalgene® jar. The use of leak proof glass and Nalgene® jars helps to enhance the safety of all personnel and preserve the evidence. Each of the Sirchie Liquid Sampling Kits is designed to have a glass jar that fits perfectly into a Nalgene® jar with pipettes for collection of the substance.

### 30ml Liquid Sampling Kit

This kit includes 30ml Nalgene® Bottles (12 each), 4ml Glass Vials (12 each), and Disposable Pipettes (12 each).



**No. LSK100**

**30ml Liquid Sampling Kit**

### 60ml Liquid Sampling Kit

This kit includes 60ml Nalgene® Bottles (12 each), 12ml Glass Vials (12 each), and Disposable Pipettes (12 each).



**No. LSK200**

**60ml Liquid Sampling Kit**

### 125ml Liquid Sampling Kit

This kit includes 125ml Nalgene® Bottles (12 each), 16ml Glass Vials (12 each), and Disposable Pipettes (12 each).



**No. LSK300**

**125ml Liquid Sampling Kit**

### 250ml Liquid Sampling Kit

This kit includes 250ml Nalgene® Bottles (12 each), 30ml Glass Vials (12 each), and Disposable Pipettes (12 each).



**No. LSK400**

**250ml Liquid Sampling Kit**

### Disposable Liquid Sampling Pipettes

Use these disposable pipettes for sampling unknown liquids. Store collected liquids in glass vials placed inside Nalgene bottles.



**No. PIP25**

**Disposable Liquid Sampling Pipettes, 25 ea.**

## Translucent Nalgene® Bottles

You have the option to get the wide-mouth Nalgene® bottles separately.



**No. NAJAR0001** Translucent Nalgene® Bottles, 30ml, 12 ea.

**No. NAJAR0002** Translucent Nalgene® Bottles, 60ml, 12 ea.

**No. NAJAR0004** Translucent Nalgene® Bottles, 125ml, 12 ea.

**No. NAJAR0008** Translucent Nalgene® Bottles, 250ml, 12 ea.

**No. NAJAR0016** Translucent Nalgene® Bottles, 500ml, 12 ea.

**No. NAJAR0032** Translucent Nalgene® Bottles, 1 liter, 6 ea.

## Amber Nalgene® Lab Quality Bottles

Store, ship and package light-sensitive liquid using these narrow-mouth Amber HDPE lab quality Nalgene® bottles. Featuring excellent chemical resistance to most acids, bases, and alcohols.



**No. NABOT9125** Amber Nalgene® Bottles, 4ml, 332 ea.

**No. NABOT9025** Amber Nalgene® Bottles, 8ml, 12 ea.

**No. NABOT9050** Amber Nalgene® Bottles, 15ml, 150 ea.

**No. NABOT0001** Amber Nalgene® Bottles, 30ml, 12 ea.

**No. NABOT0002** Amber Nalgene® Bottles, 60ml, 12 ea.

**No. NABOT0004** Amber Nalgene® Bottles, 125ml, 12 ea.

**No. NABOT0008** Amber Nalgene® Bottles, 250ml, 12 ea.

**No. NABOT0016** Amber Nalgene® Bottles, 500ml, 12 ea.

**No. NABOT0032** Amber Nalgene® Bottles, 1 liter, 6 ea.

**No. NABOT0010** Amber Nalgene® Bottles, 1 gal, 1 ea.



## NARKsafe™ Illicit Substance Testing Workstation

Fentanyl-related deaths are increasing as the opioid epidemic increases worldwide. Fentanyl is deadly if swallowed, inhaled, or absorbed through the skin, so being protected when handling suspected substances is paramount. NARKsafe™ fume hood enclosures are ideal for protecting professionals during the forensic analysis of fentanyl and other dangerous substances. Only this system offers a turbulence-free, horizontal laminar airflow with dual HEPA filtration.




NARKsafe™ enclosure provides a negative pressure environment to pull any particulate away from the operator and eliminate end-user exposure. It also allows for the addition of a bonded carbon filter, so analysis of other chemical solutions or marijuana exhibits can be safely handled in the same workspace. Available in three widths: 36", 48", 72".

### Why NARKsafe™ is the only enclosure to use when testing illicit substances:

- Horizontal laminar airflow pattern
- Dual HEPA filtration
- HEPA filters with 99.997% efficiency rating for collecting 0.3-micron particles
- Electrostatically charged pre-filter to collect larger particulate and minimize clogging of the HEPA filter
- ASHRAE 110 testing adjustable face velocity
- Ability to replace contaminated filters under negative pressure
- Blue flooring to see substances clearly when testing



**i** When testing illicit substances and changing filters, the operator should always follow Drug Enforcement Administration (DEA) and the Centers for Disease Control (CDC) guidelines for handling fentanyl, carfentanil and acetyl-fentanyl.









No. ACDD760	NARKsafe™ 36" Powder Safety Workstation, 110V AC w/prefilter, carbon filter & HEPA filters
No. ACDD760220	NARKsafe™ 36" Powder Safety Workstation, 220V AC w/prefilter, carbon filter & HEPA filters
No. ACDD770	NARKsafe™ 48" Powder Safety Workstation, 110V AC w/prefilter, carbon filter & HEPA filters
No. ACDD770220	NARKsafe™ 48" Powder Safety Workstation, 220V AC w/prefilter, carbon filter & HEPA filters
No. ACDD780	NARKsafe™ 72" Powder Safety Workstation, 110V AC w/prefilter, carbon filter & HEPA filters
No. ACDD780220	NARKsafe™ 72" Powder Safety Workstation, 220V AC w/prefilter, carbon filter & HEPA filters

## NARK II Presumptive Drug Testing Reagent Pouches

Our NARK drug identification pouches can presumptively identify several families of substances that could be abused drugs. Each of the tests is comprised of one or more chemical reagents. When a color or series of colors occur within a specific testing sequence, a positive identification may be presumed.



### Most popular reagent tests:

Catalog No.	Description	Presumptively tests for
NARK2005	Duquenois-Levine Reagent, 10/box	Marijuana, Hash, Hash Oil, THC Concentrates 
NARK2007	Scott Reagent (Modified), 10/box	Cocaine Salts (HCl), Cocaine Base (Crack) 
NARK20011	Mecke's Modified Reagent, 10/box	Heroin (White, Brown, Black Tar) 
NARK20015	Methamphetamine/MDMA Reagent, 10/box	Methamphetamine, MDMA (Ecstasy) 
NARK20023	Synthetic Cannabinoids Reagent, 5/box	Synthetic Cannabinoids (indole formulations) 
NARK20033	Fentanyl Reagent, 10/box	Fentanyl, Acetyl-fentanyl, Heroin 

Catalog No.	Description	Presumptively tests for
NARK2001	Marquis Reagent, 10/box	Opium Alkaloids (Heroin/Morphine), MDMA (Ecstasy), Amphetamines/Meth
NARK2002	Nitric Acid Reagent, 10/box	Differentiate Heroin/Morphine
NARK2003	Dille-Koppanyi Reagent, 10/box	Barbiturates
NARK2004	Ehrlich's Reagent, 10/box	LSD
NARK2008	Methadone Reagent, 10/box	Methadone
NARK2009	PCP/Methaqualone Reagent, 10/box	PCP, Methaqualone
NARK20010	Special Opiates Reagent, 10/box	Codeine, Heroin, Morphine, Oxycodone, Buprenorphine
NARK20012	Froehde's Reagent, 10/box	Talwin® (Pentazocine), DMT (Dimethyltryptamine)
NARK20013	Ephedrine Reagent, 10/box	Ephedrine, Pseudoephedrine
NARK20014	Valium Reagent, 10/box	Valium, Rohypnol, Ketamine
NARK20019	Mayer's Reagent, 10/box	General Screening, Narcotic Alkaloids
NARK20020	KN Reagent (Fast Blue B Salts), 10/box	Marijuana, Hashish, Hash Oil, THC Concentrates (DAB, WAX, BHO)
NARK20021	GHB Reagent, 5/box	GHB (Gamma Hydroxybutyrate)
NARK20022	Mandelin Reagent, 10/box	General screening for Methadone, Amphetamine-Type Compounds
NARK20024	MDPV Reagent, 10/box	MDPV, Butylone, Ethylone (bk-MDEA), Methylone, Naphyrone, Pentylone
NARK20025	Mephedrone Reagent, 10/box	Mephedrone (Synthetic Cathinones)
NARK20026	a-PVP Reagent, 10/box	a-PVP (Flakka)
NARK20029	2C Reagent, 10/box	2C-B, 2C-B-NBOMe, 2C-C, 2C-C-NBOMe, 2C-E, 2C-E-NBOMe, 2C-I, 2C-I-NBOMe, 2C-N, 2C-T-7
NARK20030	Psilocybin/Psilocin Reagent, 10/box	Psilocybin, Psilocin (Mushrooms)
NARK20031	Liebermann Reagent, 10/box	Methamphetamine, Morphine
NARK20032	Mollies Reagent, 10/box	Mollies

Catalog No.	Description
NARK2006	Acid Neutralizer Liquid, 1 oz bottle
NARK200N	Acid Neutralizer Powder, 0.5 oz. pouch, 10 each



**Pricing information:**  
For current pricing visit  
[www.sirchie.com](http://www.sirchie.com), or  
call (800) 356-7311



## NARK II Notebook and Patrol Kits

Carry the most commonly used tests with you in a protected case, keeping the glass ampoules in the pouches safe. We designed our Patrol Kit with the reagents most often encountered by patrol officers.



**No. NARK12PAK**

**NARK II 12-Pak Notebook Kit**

**No. NARK24PAK**

**NARK II 24-Pak Notebook Kit**

**No. NARK200VC**

**NARK II Patrol Kit**

## NARK II Presumptive Drug Test Kits

Popular tests packaged in convenient kits. We offer a variety of combination field tests in rugged plastic cases from 16 to 130 individual tests. Custom kits are available.



**No. NARK200CD**

**NARK II Customs/Douanes Kit**

**No. NARK200MEG**

**NARK II Meg Kit**

**No. NARK200PK**

**NARK II Port-A-Kit**

## Mollies Drug Reagent Kit

The Mollies Kit contains 50 individual tests to help you presumptively identify substances commonly known as Mollies on the street.



**No. NARK200MLE**

**NARK II Mollies Drug Reagent Kit**

## Master Synthetics Drug Reagent Kit

With this kit, you'll have all our reagent kits to conduct a presumptive drug identification for synthetic substances in one easy and convenient package. There's an acid neutralizer included so presumptive tests can be neutralized and you can dispose of the test according to your department protocols.



**No. NARK200SYN**

**NARK II Master Synthetic Drug Reagent Kit**

## NARK II Corrections Kit

Ideal for Corrections Officers, we developed this kit to include the presumptive drug analysis tests for substances commonly abused in prisons. To make it easier for you, we've included a chart to properly identify drugs during tests.



**No. NARK200CR**

**NARK II Corrections Kit**

## NARK II Cocaine ID Swipes

Presumptively identify cocaine with only trace amounts left on surfaces. Easy-to-use and results are immediate: the swipe will turn blue if the results are presumptively positive.



**No. NS20150**

**NARK II Cocaine ID Swipes**

## Thin Layer Chromatography of Marijuana Kit

No false positives will be encountered. This is a rapid, simplified, specific test for marijuana using a laboratory technique that any officer may perform.



**No. MZ020**

**Thin Layer Chromatography of Marijuana Kit**

**No. MZ020A**

**Replacement Test, 1 ea.**

## Weighing Scales

Use these scales for precise weighing of evidence. Choose from our Basic Compact Scale, Advanced Scale or Precision Balance.



**No. SCL200**

**Basic Compact Scale, +/- 0.1, 110V**

**No. SCL5000**

**Basic Compact Scale, +/- 2g, 110V**

**No. SCL5000E**

**Basic Compact Scale, +/- 2g, 220V**

**No. SCL2600**

**Advanced Scale, 2600g +/- 2g, 110V**

**No. SCL1000**

**Precision Balance, 1000g +/-0.01g, 110V**

**No. SCL1000E**

**Precision Balance, 1000g +/-0.01g, 220V**



**Often Bought Together:**  
Indicates items often  
purchased in conjunction with  
other items.




## NARK Test Reagents and Components

These reagents are in disposable tube format with each reagent in a glass ampoule.

### Our top 3 drug reagents:

- NARK Meckes Modified Reagent (Heroin.) No. NAR10016
- NARK Special Opiates (Codeine, Heroin, Oxycodone). No. NAR10022
- NARK MDPV Reagent (Synthetic Cathinones). No. NAR10024



No.	Description	Presumptively tests for
NAR10001	Mayer's Reagent, 10/box	Narcotic Alkaloids
NAR10002	Marquis Reagent, 10/box	Opium Alkaloids
NAR10003	Nitric Acid Reagent, 10/box	Differentiate Heroin/Morphine
NAR10004	Cobalt-Thiocyanate Reagent, 10/box	Cocaine/Crack
NAR10005	Dille-Koppanyi Reagent, 10/box	Barbiturates
NAR10006	Mandelin Reagent, 10/box	Amphetamines
NAR10007	Ehrlich's Reagent, 10/box	LSD
NAR10008	Duquenois Reagent, 10/box	Marijuana, Hash, Hash Oil, THC Concentrates
NAR10009	KN Reagent, 10/box	Marijuana (green plant & seeds)
NAR10013	Modified Cobalt-Thiocyanate Reagent, 10/box	Cocaine Freebase, Crack
NAR10014	Valium Reagent, 10/box	Valium, Rohypnol, Methcathinone, Ketamine
NAR10015	Meth/MDMA Reagent, 10/box	Methamphetamine, MDMA (Ecstasy)
 NAR10016	Mecke's Modified Reagent, 10/box	Heroin (white, brown, black tar)
NAR10018	Froehde's Reagent, 10/box	Talwin® (Pentazocine), DMT (Dimethyltryptamine)
NAR10019	Ephedrine Reagent, 10/box	Ephedrine, Pseudoephedrine
 NAR10022	Special Opiates Reagent, 10/box	Codeine, Heroin, Oxycodone
 NAR10024	MDPV Reagent, 10/box	MDPV (Synthetic Cathinones)
NAR10025	Mephedrone Reagent, 10/box	Mephedrone (Synthetic Cathinones)
NAR10028	GHB Reagent, 5/box	GHB (Gamma Hydroxybutyrate)
NAR10029	PCP/Methaqualone Reagent, 10/box	PCP, Methaqualone
NAR10030	Psilocybin/Psilocin Reagent, 10/box	Psilocybin, Psilocin (mushrooms)
NAR10010	Polyester Bags, 10/box	
NAR10011	Acid Neutralizer, 1 oz. Bottle	

## NARK Analysis Reagent Kit

We included reagents for the most frequently abused drugs in this field test kit so you always have them on hand. Each individual reagent is sealed in a thin, glass ampoule which eliminates cross-contamination with other reagents.

Reagents included for: Marijuana, GHB, Hashish, THC, Cocaine, Crack/Freebase, Amphetamines, Barbiturates, Narcotic Alkaloids, Opium Alkaloids, PCP, Hallucinogens, Heroin, Morphine.

**No. NAR100**

**NARK Analysis Reagent Kit**

