

Hotplates, Stirrers, Rockers, and Scales

Laboratory Equipment

HOTPLATES, STIRRERS, AND HOTPLATE/STIRRERS

These hotplates, stirrers and hotplate-stirrers feature an exceptionally durable, chemical resistant white ceramic work surface. Their space-efficient design (8 X 9 in. footprint), makes them ideal for use on crowded benchtops and inside of bio-hoods.

Advanced microprocessor controls with convenient turn knobs allow quick, precise adjustment and maintenance of speed and temperature. Safety indicator LED's on the front panel indicate when the heating and/or stirring functions have been activated. With a square 7.5 in work surface, all three models are compatible with a wide variety of popular sizes of borosilicate glass beakers, flasks, bottles and other vessels. The hotplate/magnetic stirrer model (H4000-HS) comes complete with a support rod for mounting thermometers and temperature probes.

BENCH ROCKER

The powerful motor and heavy-duty drive system remain reliable and consistent, even under constant usage in a cold room (min. 0°C) or incubator (max. 65°C). The rocking speed and tilt angle can easily be adjusted to vary between gentle and aggressive agitation. Both the standard platform and the (optional) stacking platform are constructed of stainless steel, preventing against any rusting or corrosion.

CATALOG NO.	DESCRIPTION	PRICE
HP110	Hotplate, 110V	
HP220	Hotplate, 220V	
SP110	Magnetic Stirrer, 110V	
SP220	Magnetic Stirrer, 220V	
HPS110	Hotplate / Stirrer Combination, 110V	
HPS220	Hotplate / Stirrer Combination, 220V	
BR110	Bench Rocker, 16" x 12", 110V	
BR220	Bench Rocker 16" x 12", 220V	

BASIC COMPACT SCALES

The Basic Compact Scales offer all of the necessary features for quick and easy weighing. Lightweight and portable these scales feature battery as well as AC operation, are splashproof and offer 3 weighing units (g, lb, oz).

ADVANCED SCALES

The Advanced Scale is the best choice for simple operation and economy, featuring a stainless steel pan, three-point overload protection and battery as well as AC power operation.

PRECISION BALANCES

Precision balances have what it takes for basic lab work, field use and various industrial applications. These balances feature 15 weighing units, three-point overload protection, RS-232 and USB interface, internal and external calibration and battery as well as AC operation.

CATALOG NO.	DESCRIPTION	PRICE
SCL200	Basic Compact Scale, 200g +/- 0.1, 110V	
SCL200E	Basic Compact Scale, 200g +/- 0.1, 220V	
SCL5000	Basic Compact Scale, 5000g +/- 2g, 110V	
SCL5000E	Basic Compact Scale, 5000g +/- 2g, 220V	
SCL2600	Advanced Scale, 2600g +/- 0.1g, 110V	
SCL2600E	Advanced Scale, 2600g +/- 0.1g, 220V	
SCL1000	Precision Balance, 1000g +/- 0.01g, 110V	
SCL1000E	Precision Balance, 1000g +/- 0.01g, 220V	

NEW!



Hotplate

NEW!

Magnetic Stirrer



NEW!



Hotplate/Stirrer

NEW!



Bench Rocker

NEW!



Basic Compact Scale

NEW!

Advanced Scale



NEW!



Precision Balance

Helmets & Body Shields



Premier Crown brand helmets are assembled at Sirchie's manufacturing facility in Youngsville, NC.



Protective Body Shields

For more information on Premier Crown products call us at 800.356.7311 or 919.554.2244 or go to www.premiercrown.com

PREMIER CROWN Helmets and Body Shields

Established in 1979, Premier Crown is honored to continue providing high-quality tactical and safety equipment needed to protect the well being of our nation's law enforcement and correctional officers, as well as all law enforcement professionals around the world. We pride ourselves on providing quality products through a company-wide commitment that incorporates progressive production techniques, state-of-the-art machinery and continuous product research, development and improvement. At Premier Crown, we've developed new products that will greatly increase the capability of tactical operations for riot control squads and correctional cell extraction teams.

RIOT DUTY HELMETS



GENERAL DUTY HELMETS

