



TECHNICAL INFORMATION

MAGNUCLEI™ Magnetic Latent Print Powder

Catalog Nos. MAG100, MAG10016

INTRODUCTION

MAGNUCLEI™ is the result of years of research and development to fill the need for a dual-purpose powder for finer dusting applications. MAGNUCLEI™ powders have been milled to average particle sizes less than half the size of regular magnetic powders. It very effective on glass, plastic, wood, ceramic, and metallic surfaces—use on paper surfaces without leaving dark backgrounds. It works well on imitation and genuine leather if developed within a short time of being handled. (This powder DOES NOT interfere with Ninhydrin, DFO, or other chemical treatments.)

CAUTION

- Before using this product, consult the appropriate Material Safety Data Sheets (MSDS) found on our website at www.sirchie.com/support.



100 HUNTER PLACE, YOUNGSVILLE, NC 27596 USA

Ph: (919) 554-2244, (800) 356-7311 • Fax: (919) 554-2266, (800) 899-8181 • Web: www.sirchie.com • Email: sirchieinfo@sirchie.com

Leading Particulars

- Finely milled to a particle size of less than 20 microns (half the particle size of regular magnetic powders).
- Very effective on glass, plastic, wood, ceramic, metal surfaces, imitation or genuine leather, and paper surfaces.
- Reveal latent prints with a single pass of the magnetic wand.
- Allows for normal lifting of developed latent prints with lifting tape and hinge lifters.

NOTE: We recommend the use of Rubber-Jell™ Lifting Tape or Rubber/Gel Lifters if developed latent prints are found on rough, textured, irregular, or curved surfaces.



Lifts of resultant prints are shown on both white and black backing material.

