

SECTION 1: Identification

Identification

Product form : Substance

Substance name : LV502 Crystal Violet (Methyl Violet) Powder

Product code LV502; LV5021

Recommended use and restrictions on use

Use of the substance/mixture : Dyestuff/pigment

Supplier

SIRCHIE

100 Hunter Place

Youngsville, NC 27596 - USA

T 919-554-2244; 800-356-7311 - F 919-554-2266; 800-899-8181

http://www.sirchie.com

Emergency telephone number

: 1.800.424.9300 (USA) +1-703-527-3887 (INTL) Emergency number

CHEMTREC: 1.800.424.9300

SECTION 2: Hazard(s) identification

Classification of the substance or mixture

GHS US classification

Acute toxicity (oral) Category 4 H302 Harmful if swallowed Serious eye damage/eye irritation Category 2 H319 Causes serious eye irritation

Full text of H statements: see section 16

GHS Label elements, including precautionary statements

GHS US labeling

Hazard pictograms (GHS US)



Signal word (GHS US) : Warning

H302 - Harmful if swallowed Hazard statements (GHS US)

H319 - Causes serious eye irritation

Precautionary statements (GHS US) P264 - Wash hands, forearms and face thoroughly after handling.

P270 - Do not eat, drink or smoke when using this product.

P280 - Wear protective gloves/protective clothing/eye protection/face protection. P301+P312 - If swallowed: Call a poison center or doctor if you feel unwell.

P305+P351+P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove

contact lenses, if present and easy to do. Continue rinsing.

P330 - Rinse mouth.

P337+P313 - If eve irritation persists: Get medical advice/attention.

P501 - Dispose of contents/container to hazardous or special waste collection point, in

accordance with local, regional, national and/or international regulation.

Other hazards which do not result in classification

Other hazards not contributing to the : None under normal conditions.

classification

2.4. **Unknown acute toxicity (GHS US)**

Not applicable

SECTION 3: Composition/Information on ingredients

3.1. **Substances**

: LV502 Crystal Violet (Methyl Violet) Powder Name

02/16/2021 EN (English US) Page 1

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Name	Product identifier	%	GHS US classification
C.I. basic violet 1	(CAS-No.) 8004-87-3	75	Acute Tox. 4 (Oral), H302
			Eye Irrit. 2A, H319

Full text of hazard classes and H-statements : see section 16

Not applicable

SECTION 4: First-aid measures

Description of first aid measures

: Call a poison center/doctor/physician if you feel unwell. First-aid measures general

: Remove the victim into fresh air. If you feel unwell, seek medical advice. First-aid measures after inhalation

First-aid measures after skin contact Wash immediately with lots of water (15 minutes)/shower. Take victim to a doctor/medical

service if irritation persists.

: Rinse immediately with plenty of water for 15 minutes. Get medical advice/attention. First-aid measures after eye contact

First-aid measures after ingestion Do not induce vomiting. Immediately call a poison center or doctor/physician.

Most important symptoms and effects (acute and delayed)

Potential Adverse human health effects and

: Based on available data, the classification criteria are not met.

symptoms

Symptoms/effects after skin contact : May cause moderate irritation. Symptoms/effects after eye contact : Causes serious eye damage. Symptoms/effects after ingestion : May be harmful if swallowed.

Immediate medical attention and special treatment, if necessary

No additional information available

SECTION 5: Fire-fighting measures

Suitable (and unsuitable) extinguishing media

Suitable extinguishing media : Adapt extinguishing media to the environment for surrounding fires.

Specific hazards arising from the chemical

Reactivity in case of fire : No data available.

Special protective equipment and precautions for fire-fighters

Firefighting instructions : Use water spray or fog for cooling exposed containers. Exercise caution when fighting any

chemical fire. Prevent fire-fighting water from entering environment.

Protection during firefighting : Do not enter fire area without proper protective equipment, including respiratory protection.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Emergency procedures : Evacuate unnecessary personnel.

6.1.2. For emergency responders

Protective equipment : Equip cleanup crew with proper protection.

Emergency procedures Ventilate area.

Environmental precautions

Prevent entry to sewers and public waters. Notify authorities if liquid enters sewers or public waters.

Methods and material for containment and cleaning up

: On land, sweep or shovel into suitable containers. Minimize generation of dust. Store away Methods for cleaning up

from other materials.

Reference to other sections

See Heading 8. Exposure controls and personal protection.

SECTION 7: Handling and storage

Precautions for safe handling

Precautions for safe handling : Avoid dust formation. Avoid contact with skin and eyes. Do not breathe dust.

02/16/2021 EN (English US) 2/6

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Keep container closed when not in use. Store in original container.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

LV502 Crystal Violet (Methyl Violet) Powder

No additional information available

C.I. basic violet 1 (8004-87-3)

No additional information available

8.2. Appropriate engineering controls

8.3. Individual protection measures/Personal protective equipment

Personal protective equipment:

Dust formation: dust mask. Protective goggles. Gloves.

Hand protection:

Wear protective gloves.

Eye protection:

Chemical goggles or safety glasses

Respiratory protection:

Wear appropriate mask

Personal protective equipment symbol(s):







Other information:

Do not eat, drink or smoke during use.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Solid
Appearance : Powders.
Color : Blue-green
Odor : odorless

Odor threshold : No data available : No data available рΗ : No data available Melting point : No data available Freezing point : No data available Boiling point Flash point : No data available Relative evaporation rate (butyl acetate=1) : No data available : Non flammable. Flammability (solid, gas) Vapor pressure No data available : No data available Relative vapor density at 20 °C : No data available Relative density

02/16/2021 EN (English US) 3/6

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

: Soluble in water. Solubility Partition coefficient n-octanol/water (Log Pow) : No data available Auto-ignition temperature No data available Decomposition temperature : No data available Viscosity, kinematic : No data available Viscosity, dynamic : No data available **Explosion limits** : No data available Explosive properties : No data available Oxidizing properties : No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

No data available.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

Not established.

10.4. Conditions to avoid

Direct sunlight. Extremely high or low temperatures.

10.5. Incompatible materials

Strong acids. Strong bases.

10.6. Hazardous decomposition products

fume. Carbon monoxide. Carbon dioxide.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity (oral) : Harmful if swallowed.

Acute toxicity (dermal) : Not classified

Acute toxicity (inhalation) : Not classified

ATE US (oral) 550.667 mg/kg body weight

C.I. basic violet 1 (8004-87-3)

LD50 oral rat 413 mg/kg (Rat, Oral)

Skin corrosion/irritation : Not classified

Serious eye damage/irritation : Causes serious eye irritation.

Respiratory or skin sensitization : Not classified
Germ cell mutagenicity : Not classified
Carcinogenicity : Not classified
Reproductive toxicity : Not classified

Specific target organ toxicity – single exposure : Not classified

Specific target organ toxicity - repeated : Not classified

exposure

Aspiration hazard : Not classified
Viscosity, kinematic : No data available

Potential Adverse human health effects and

symptoms

: Based on available data, the classification criteria are not met.

Symptoms/effects after skin contact : May cause moderate irritation. Symptoms/effects after eye contact : Causes serious eye damage.

02/16/2021 EN (English US) 4/6

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Symptoms/effects after ingestion : May be harmful if swallowed.

SECTION 12: Ecological information

12.1. Toxicity

No additional information available

12.2. Persistence and degradability

LV502 Crystal Violet (Methyl Violet) Powder		
Persistence and degradability Not established.		
C.I. basic violet 1 (8004-87-3)		
Persistence and degradability	Biodegradability in water: no data available.	

12.3. Bioaccumulative potential

LV502 Crystal Violet (Methyl Violet) Powder		
Bioaccumulative potential	Not established.	
C.I. basic violet 1 (8004-87-3)		
Bioaccumulative potential	No bioaccumulation data available.	

12.4. Mobility in soil

No additional information available

12.5. Other adverse effects

Other information : Avoid release to the environment.

SECTION 13: Disposal considerations

13.1. Disposal methods

Product/Packaging disposal recommendations : Dispose in a safe manner in accordance with local/national regulations.

Ecology - waste materials : Avoid release to the environment.

SECTION 14: Transport information

Department of Transportation (DOT)

In accordance with DOT

Other information : No supplementary information available.

Transportation of Dangerous Goods

Transport by sea

Air transport

SECTION 15: Regulatory information

15.1. US Federal regulations

No additional information available

15.2. International regulations

CANADA

No additional information available

EU-Regulations

No additional information available

National regulations

No additional information available

15.3. US State regulations

No additional information available

SECTION 16: Other information

02/16/2021 EN (English US) 5/6

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

: REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE Data sources

COUNCIL of 16 December 2008 on classification, labeling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending

Regulation (EC) No 1907/2006.

Training advice : Normal use of this product shall imply use in accordance with the instructions on the packaging.

Other information None.

Full text of H-phrases:

	H302	Harmful if swallowed	
	H319	Causes serious eye irritation	
NFF	PA health hazard	: 1 - Materials that, under emergency conditions, can cause significant irritation.	
NFF	PA fire hazard	: 2 - Materials that must be moderately heated or exposed to relatively high ambient temperatures before ignition can occur.	
NFF	PA reactivity	: 0 - Material that in themselves are normally stable, even	



Health : 1 Slight Hazard - Irritation or minor reversible injury possible

under fire conditions.

1 Slight Hazard - Materials that must be preheated before ignition will occur. Includes liquids, Flammability

solids and semi solids having a flash point above 200 F. (Class IIIB)

Physical : 0 Minimal Hazard - Materials that are normally stable, even under fire conditions, and will NOT

react with water, polymerize, decompose, condense, or self-react. Non-Explosives.

Personal protection

E - Safety glasses, Gloves, Dust respirator

SDS US (GHS HazCom 2012)

The information above is believed to be accurate and represents the best information currently available to us. However, we make no warranty of merchantability or any other warranty, expressed or implied, with respect to such information, and we assume no liability resulting from its use. Users should make their own investigation to determine the suitability of the information for their particular purposes.

02/16/2021 EN (English US) 6/6